Application No. 10/073,899 Attorney Docket No. 07414.0054-00

 (Amended) A data collection method for scanning a scan window comprising one or more channels comprising the steps of:

detecting an integrated signal (S) across [a] the scan window comprising one or more channels using an integrating detector; and

calculating a velocity-normalized integrated signal (Sn) as a function of a scan velocity and the integrated signal S.

3. (Amended) The method of claim 1 wherein the step of calculating the velocity-normalized integrated signal (Sn) comprises:

measuring a channel width (w);

determining a time for traversing the channel width (t); and

computing [a] the velocity-normalized integrated signal according to the equation Sn=S/(w/t).

- 4. (Amended) The method of claim 1 wherein the step of calculating the velocity-normalized integrated signal (Sn) comprises subtracting a detector offset So from [an] the integrated signal (S).
- 16. (Amended) An apparatus for scanning a plurality of channels comprising: means for detecting an integrated signal (S) across a scan window comprising [one or more] the plurality of channels using an integrating detector; and

computer means for receiving the integrated signal S and determining a scan velocity and for calculating a velocity-normalized integrated signal (Sn) as a function of

FINNEGAN HENDERSON FARABOW GARRETT & DUNNER LLP

1300 I Street, NW Washington, DC 20005 202.408.4000 Fax 202.408.4400 www.finnegan.com

Application No. 10/073,899 Attorney Docket No. 07414.0054-00

the scan velocity and the integrated signal S.

17. (Twice Amended) An apparatus for scanning a scan window having one or more channels comprising:

an integrating detector;

a scanner for effecting a scanning of the integrating detector relative to [a] the scan window comprising one or more channels, wherein an integrated signal (S) is detected by scanning the integrating detector relative to the scan window; and

a computer for receiving the integrated signal S and for determining a scan velocity and for calculating a velocity-normalized integrated signal (Sn).

23. (Amended) The program storage device of claim 21 wherein the step of calculating the velocity-normalized integrated signal (Sn) comprises:

measuring a channel width (w);

determining a time for traversing the channel width (t); and computing [a] the velocity-normalized integrated signal according to the equation Sn=S/(w/t).

- 25. (Amended) The program storage device of claim 24 wherein [a] the channel width (w) is measured by counting steps in the stepper motor.
- 31. (Amended) A data collection method for scanning a scan window comprising:

FINNEGAN HENDERSON FARABOW GARRETT & DUNNERL

1300 I Street, NW Washington, DC 20005 202.408.4000 Fax 202.408.4400 www.finnegan.com

Application No. 10/073,899 Attorney Docket No. 07414.0054-00

detecting an integrated signal (S) across [a] the scan window comprising one or more channels using an integrating detector;

determining an integration time (ti) for the integrated signal; and calculating a velocity-normalized integrated signal (Sn), the calculating comprising dividing the integrated signal (S) by the integration time (ti).

- 45. (Amended) The apparatus of claim 16, further comprising means for determining an integration time (ti) for the integrated signal; and wherein the calculating the velocity-normalized <u>integrated</u> signal comprises dividing the integrated signal (S) by the integration time (ti).
- 52. (Amended) The apparatus of claim 17, further comprising a timer configured to determine an integration time (ti) for the integrated signal; and wherein the calculating the velocity-normalized <u>integrated</u> signal comprises dividing the integrated signal (S) by the integration time (ti), and the scan window comprises more than one channel.

## <u>REMARKS</u>

Applicant acknowledges the Examiner's indication that claims 28-65 contain allowable subject matter. (Office Action, page 8)

FINNEGAN HENDERSON FARABOW GARRETT & DUNNER LLP

1300 I Street, NW Washington, DC 20005 202.408.4000 Fax 202.408.4400 www.finnegan.com